



SOL-MILLENNIUM MEDICAL

Product Catalog

www.sol-m.com

Global v1.3

SOL-MILLENNIUM®
Building a healthier tomorrow

CONTENT



COMPANY INTRODUCTION	3
NEW & AVAILABLE	4
General Use	7
SOL-GUARD™ Safety Syringe	8
SOL-CARE™ Safety Needle	9
SOL-CARE™ Safety Needle and Syringe Combinations	9
SOL-CARE™ Safety Syringe with Fixed Needle	11
SOL-CARE™ Luer Lock Safety Syringe with Exchangeable Needle	12
SOL-CARE™ Luer Lock Safety Syringe without needle	12
SOL-M™ Traditional Syringe with and without needle	14
SOL-M™ Hypodermic Needle	16
Diabetes	17
SOL-GUARD™ Insulin Syringe	18
SOL-CARE™ Insulin Safety Syringe with Fixed Needle	19
SOL-M™ Insulin Syringe with Fixed Needle	20
SOL-M™ Pen Needle	21
Allergy	22
SOL-CARE™ Safety Syringe Allergy Tray	23
Immunization	24
SOL-CARE™ Safety Syringe Flu Tray	25
SOL-CARE™ Safety Syringe for Immunization	26
SOL-Millennium® Safety Box	27
Pharmacy	28
SOL-M™ Blunt Fill Needle	29
SOL-M™ Blunt Fill Needle with 5 micron filter	29
Safety Products for Nuclear Medicine	30
SOL-CARE™ Luer Lock Safety Syringe	31
- Fits most Nuclear Medicine Shields	31
Blood Collection	32
SOL-CARE™ Safety Blood Collection Tube Holder Set	33
Oral Feeding	34
SOL-M™ Oral Dispensing Syringe	35
Vetenerian	36
SOL-M™ Insulin Syringe with Fixed Needle For U40 and U100	37

COMPANY INTRODUCTION

Dedicated to Building a Healthier Tomorrow, improve the safety, and well-being of healthcare workers and their patients.

Sol-Millennium Group is a Global Healthcare Company providing innovative products and value added services through product development, manufacturing and distribution focusing on safe injection and disposal devices.

Marketing and Supplying the global market through five operating entities:

- * Sol-Millennium Medical Products Co., Ltd.
- * Sol-Millennium Aust Pty.
- * Sol-Millennium Medical, Inc.
- * Sol-Millennium Brasil Imp.Exp. S/A
- * Sol-Millennium Europe Sp. z o.o.

SOL-MILLENNIUM

Building a healthier tomorrow

A New Choice for Safety Injection by **SOL-CARE™ Safety Needle**

SOL-CARE™ Safety Needle is compatible with all luer lock and luer slip standard syringes, easy to use, protecting the healthcare worker from potential in needlestick injury.



*Single-handed activation
Needle completely covered*

.....

*For detail product information,
go to page 9.*

**NOW
AVAILABLE**

SOL-CARE™

Safety Blood Collection

Tube Holder Set

**NEW
AVAILABLE**

*Integrated safety mechanism
compatible with standard multi-sample needles
.....*

*For detail product information,
go to page 20.*




Complies with OSHA regulations to help eliminate cross contamination.

SOL-GUARD™ Safety Syringe

is a safety-engineered syringe for various injections. It is designed to reduce the risk of needlestick injuries in medication transportation and injection process for healthcare workers.

NEW COMING SOON



*Transparent Locking Safety Shield
Transport Position
Low Dead Space
.....*

*For detail product information,
go to page 7 & 16 .*

NEW COMING SOON

SOL-M™ Insulin Syringe with Fixed Needle For U40 Insulin

Sol-Millennium offers innovative products for the Vet Market for comfort and convenience for diabetes care for your pets.

NEW COMING SOON



New Veterinerian Syringe

*Low dead space
Optimal needle sharpness
Transparent barrel with bold scale
.....*

*For detail product information,
go to page 8.*



/ **General Use**

SOL-GUARD™


Safety Syringe



SOL-GUARD™ Safety Syringe, engineered safety mechanism to potentially reduce the risk of needlestick injury. The SOL-GUARD™ Safety Syringe, has a transport mode, clear graduations for visibility of draw-up of medication and low dead space.

Features & Benefits

- **Special Design Safety Shield** - Protection to healthcare worker during transporting device, and after use. Easy to lock either right or left twist, confirmation for final locking position is an audible click
- **Transport Position** - Potentially reducing the risk of needle stick injury while moving locations
- **Low Dead Space** - Minimizes medication waste
- **Clear, Bold Graduations** - Accuracy in draw up of medication



REF	Description	Qty
-----	-------------	-----

1mL SOL-GUARD™ Safety Syringe

BLISTER PACKAGE

2102215	1ml	22G x 1 1/2"	100/Box; 800/Case
200034SG	1ml	25G x 1"	100/Box; 800/Case
200018SG	1ml	25G x 5/8"	100/Box; 800/Case
200019SG	1ml	27G x 1/2"	100/Box; 800/Case
200067SG	1ml	28G x 1/2"	100/Box; 800/Case

1.5mL SOL-GUARD™ Safety Syringe

BLISTER PACKAGE

2152215	1.5ml	22G x 1 1/2"	100/Box; 800/Case
---------	-------	--------------	-------------------

Please contact us for the availability in your country.

SOL-CARE™

Safety Needle

Safety Needle and Syringe Combinations




Features & Benefits

- **Easy to use** - one-hand activation, minimal training required
- **Safety mechanism is an integrated part of the needle** - needle remains secure in double locking mechanism
- **Bevel-up needle** - Allows for low-angle injections
- **Easy positioned safety-arm** - Optimum visualization of injection site
- **Patent Locking design** - Needle remains securely covered
- **Needle hub** - compatible with all standard luer lock and luer slip tip syringes
- **Low dead space** - Ability to minimize the medication waste
- **Audible "click"** - Confirmation safety mechanism has been activated
- **ISO color coded** - Easy and quick identification of the needle gauge
- **Complete Product Selection** - Includes a wide range of sizes to meet clinical requirements

Two Options to activate:



OR



REF	Description	Qty
SOL-CARE™ Safety Needles		
SN3012	30 G x 1/2 in.	100/Inner Box; 800/Case
SN2758	27 G x 5/8 in.	100/Inner Box; 800/Case
SN2712	27 G x 1/2 in.	100/Inner Box; 800/Case
SN2734	27 G x 3/4 in.	100/Inner Box; 800/Case
SN2715	27 G x 1 1/2 in.	100/Inner Box; 800/Case
SN27125	27 G x 1 1/4 in.	100/Inner Box; 800/Case
SN2638	26 G x 3/8 in.	100/Inner Box; 800/Case
SN2612	26 G x 1/2 in.	100/Inner Box; 800/Case
SN2558	25 G x 5/8 in.	100/Inner Box; 800/Case
SN2510	25 G x 1 in.	100/Inner Box; 800/Case
SN2515	25 G x 1 1/2 in.	100/Inner Box; 800/Case
SN2358	23 G x 5/8 in.	100/Inner Box; 800/Case
SN2310	23 G x 1 in.	100/Inner Box; 800/Case
SN23125	23 G x 1 1/4 in.	100/Inner Box; 800/Case
SN2315	23 G x 1 1/2 in.	100/Inner Box; 800/Case
SN2210	22 G x 1 in.	100/Inner Box; 800/Case
SN22125	22 G x 1 1/4 in.	100/Inner Box; 800/Case
SN2215	22 G x 1 1/2 in.	100/Inner Box; 800/Case
SN2158	21 G x 5/8 in.	100/Inner Box; 800/Case
SN2110	21 G x 1 in.	100/Inner Box; 800/Case
SN21125	21 G x 1 1/4 in.	100/Inner Box; 800/Case
SN2115	21 G x 1 1/2 in.	100/Inner Box; 800/Case
SN2010	20 G x 1 in.	100/Inner Box; 800/Case
SN2015	20 G x 1 1/2 in.	100/Inner Box; 800/Case
SN1910	19 G x 1 in.	100/Inner Box; 800/Case
SN1915	19 G x 1 1/2 in.	100/Inner Box; 800/Case
SN1810	18 G x 1 in.	100/Inner Box; 800/Case
SN1815	18 G x 1 1/2 in.	100/Inner Box; 800/Case
SN1820*	18 G x 50mm	100/Inner Box; 800/Case
SN1920*	19 G x 50mm	100/Inner Box; 800/Case
SN2005*	20 G x 50mm	100/Inner Box; 800/Case
SN2105*	21 G x 50mm	100/Inner Box; 800/Case
SN2205*	22 G x 50mm	100/Inner Box; 800/Case

* Not available in US

1ml SOL-CARE™ Luer Lock Syringe with Safety Needle

BLISTER PACKAGE

13005SN	1ml with 30 G x 1/2 in.	100/Inner Box; 1200/Case
12705SN	1ml with 27 G x 1/2 in.	100/Inner Box; 1200/Case
12558SN	1ml with 25 G x 5/8 in.	100/Inner Box; 1200/Case
12510SN	1ml with 25 G x 1 in.	100/Inner Box; 1200/Case

3ml SOL-CARE™ Luer Lock Syringe with Safety Needle

BLISTER PACKAGE

32558SN	3ml with 25 G x 5/8 in.	100/Inner Box; 1200/Case
32510SN	3ml with 25 G x 1 in.	100/Inner Box; 1200/Case
32310SN	3ml with 23 G x 1 in.	100/Inner Box; 1200/Case
323125SN	3ml with 23 G x 1 1/4 in.	100/Inner Box; 1200/Case
32315SN	3ml with 23 G x 1 1/2 in.	100/Inner Box; 1200/Case
32210SN	3ml with 22 G x 1 in.	100/Inner Box; 1200/Case
32215SN	3ml with 22 G x 1 1/2 in.	100/Inner Box; 1200/Case
32110SN	3ml with 21 G x 1 in.	100/Inner Box; 1200/Case
32115SN	3ml with 21 G x 1 1/2 in.	100/Inner Box; 1200/Case
32010SN	3ml with 20 G x 1 in.	100/Inner Box; 1200/Case
32015SN	3ml with 20 G x 1 1/2 in.	100/Inner Box; 1200/Case

5ml SOL-CARE™ Luer Lock Syringe with Safety Needle

BLISTER PACKAGE

52215SN	5ml with 22 G x 1 1/2 in.	100/Inner Box; 1200/Case
52110SN	5ml with 21 G x 1 in.	100/Inner Box; 1200/Case
52115SN	5ml with 21 G x 1 1/2 in.	100/Inner Box; 1200/Case
52010SN	5ml with 20 G x 1 in.	100/Inner Box; 1200/Case
52015SN	5ml with 20 G x 1 1/2 in.	100/Inner Box; 1200/Case

10ml SOL-CARE™ Luer Lock Syringe with Safety Needle

BLISTER PACKAGE


102215SN	10ml with 22 G x 1 1/2 in.	50/Inner Box; 600/Case
102110SN	10ml with 21 G x 1 in.	50/Inner Box; 600/Case
102115SN	10ml with 21 G x 1 1/2 in.	50/Inner Box; 600/Case
102010SN	10ml with 20 G x 1 in.	50/Inner Box; 600/Case
102015SN	10ml with 20 G x 1 1/2 in.	50/Inner Box; 600/Case

Please contact us for the availability in your country.

SOL-CARE™

Safety Syringe with Fixed Needle

SOL-CARE™ Safety Syringes are designed to reduce the risk of needlestick injury to healthcare workers. The needle is retracted into the barrel of the syringe and secured with our patent pending locking ring eliminating syringe reuse.



Features & Benefits

- Special designed plunger with fit tip reduces dead space in syringe
- Ultra sharp needle minimizes pain on injection
- Patent Pending Locking Ring secures the needle into the barrel of the syringe

Easy to use


Fingers remain behind the needle during activation of safety mechanism

Syringe cannot be reused

Clear, bold gradations

Visual and audible confirmation of safety activation

Potential reduction in sharps disposal costs



REF	Description	Qty
-----	-------------	-----

1ml SOL-CARE™ Safety Syringe with Fixed Needle		
---	--	--

BLISTER PACKAGE

100040IM	1ml	22G x 1 1/2"	100/Box; 2400/Case
100018IM	1ml	25G x 5/8"	100/Box; 1000/Case
100034IM	1ml	25G x 1"	100/Box; 1000/Case
100067IM	1ml	28G x 1/2"	100/Box; 1000/Case
100082IM	1ml	30G x 1/2"	100/Box; 1000/Case
100019IM	1ml	27G x 1/2"	100/Box; 1000/Case
100070IM	1ml	26G x 3/8"IDB	100/Box; 1000/Case

Please contact us for the availability in your country.

SOL-CARE™

Luer Lock Safety Syringe with Exchangeable Needle Luer Lock Safety Syringe without needle

The SOL-CARE™ Safety Syringe is available with and without needles. SOL-CARE™ Safety Syringes without needle can be used any standard luer lock needle.

Features & Benefits

- All SOL-CARE™ safety syringes are now enhanced with our patented no gap design allowing clinicians to push plunger to zero line when injecting air into a vial
- Exchangeable needles allow clinician to change needle after draw-up of medication
- Special designed plunger with fit tip reduces dead space in syringe
- Patent Pending Locking Ring secures the needle in the barrel of the syringe preventing re-use of syringe

Compatible with all standard hypodermic needles

Easy to use


Fingers remain behind the needle during activation of safety mechanism










Clear, bold gradations

Visual and audible confirmation of safety activation

Potential reduction in sharps disposal costs

Ultra sharp needle minimizes injection pain












REF	Description	Qty
3ml SOL-CARE™ Luer Lock Safety Syringe with Exchangeable Needle		
BLISTER PACKAGE		
100074IM	 3ml 20G x 1"	100/Box; 800/Case
100084IM	 3ml 20G x 1 1/2"	100/Box; 800/Case
100075IM	 3ml 21G x 1"	100/Box; 800/Case
100072IM	 3ml 21G x 1 1/2"	100/Box; 800/Case
100085IM	 3ml 22G x 1"	100/Box; 800/Case
100076IM	 3ml 22G x 1 1/2"	100/Box; 800/Case
100077IM	 3ml 23G x 1"	100/Box; 800/Case
100078IM	 3ml 25G x 5/8"	100/Box; 800/Case
100079IM	 3ml 25G x 1"	100/Box; 800/Case

3ml SOL-CARE™ Luer Lock Safety Syringe without Needle

BLISTER PACKAGE			
120006IM	3ml	without Needle	100/Box; 800/Case







5ml SOL-CARE™ Luer Lock Safety Syringe with Exchangeable Needle

BLISTER PACKAGE			
140084IM	 5ml 20G x 1 1/2"	100/Box; 1200/Case	
140072IM	 5ml 21G x 1 1/2"	100/Box; 1200/Case	
140060IM	 5ml 21G x 1 1/4"	100/Box; 1200/Case	
140075IM	 5ml 21G x 1"	100/Box; 1200/Case	
140076IM	 5ml 22G x 1 1/2"	100/Box; 1200/Case	
140062IM	 5ml 22G x 1 1/4"	100/Box; 1200/Case	
140085IM	 5ml 22G x 1"	100/Box; 1200/Case	
140065IM	 5ml 23G x 1 1/4"	100/Box; 1200/Case	
140077IM	 5ml 23G x 1"	100/Box; 1200/Case	

5ml SOL-CARE™ Luer Lock Safety Syringe without Needle

BLISTER PACKAGE			
120007IM	5ml	without Needle	100/Box; 800/Case

10ml SOL-CARE™ Luer Lock Safety Syringe with Exchangeable Needle

BLISTER PACKAGE			
160084IM	 10ml 20G x 1 1/2"	100/Box; 1200/Case	
160075IM	 10ml 21G x 1"	100/Box; 1200/Case	
160072IM	 10ml 21G x 1 1/2"	100/Box; 1200/Case	
160062IM	 10ml 22G x 1 1/4"	100/Box; 1200/Case	
160085IM	 10ml 22G x 1"	100/Box; 1200/Case	
160076IM	 10ml 22G x 1 1/2"	100/Box; 1200/Case	
160065IM	 10ml 23G x 1 1/4"	100/Box; 1200/Case	
160077IM	 10ml 23G x 1"	100/Box; 1200/Case	

10ml SOL-CARE™ Luer Lock Safety Syringe without Needle

BLISTER PACKAGE			
120008IM	10ml	without Needle	100/Box; 800/Case

20ml SOL-CARE™ Luer Lock Safety Syringe without Needle

BLISTER PACKAGE			
120009IM	20ml	without Needle	50/Box; 600/Case

Please contact us for the availability in your country.

SOL-M™

Traditional Syringe

with and without needle

SOL-M™ Traditional Syringes are designed with clear barrels, bold markings, latex free and a plunger stop for prevention of plunger pull-out. SOL-M™ Traditional Syringes have low dead space to allow for accurate administration of medication and reduce medication waste. The SOL-M™ syringe can be combined with regular and intradermal needles

Features & Benefits

- **Sharp needle** - Tricut bevel minimizes pain on injection
 - **Transparent barrel and bold gradations** - Clear, precise accuracy of medication draw-up and delivery
 - **Semitransparent Needle Hub** - Clear visualization of flashback
 - **Specialized Silicon Lubrication** - Offers smooth plunger movement
- Plunger Stopper - Prevent plunger pull out from barrel

REF	Description	Qty
-----	-------------	-----

0.5ml SOL-M™ Tuberculin Syringe

BLISTER PACKAGE

180527	0.5ml 27G x 1/2"	100/Inner Box; 800/Case
180526IDB	0.5ml 26G x 3/8" IDB	100/Inner Box; 800/Case

0.5ml SOL-M™ Allergy Syringe Tray

TRAY PACKAGE

180527T	0.5ml 27G x 1/2"	25/Inner Box; 1000/Case
180527IDBT	0.5ml 27G x 3/8" IDB	25/Inner Box; 1000/Case

1ml SOL-M™ Tuberculin Syringe

BLISTER PACKAGE

181027	1ml 27G x 1/2"	100/Inner Box; 800/Case
181027S	1ml 27G x 1/2" Slip tip	100/Inner Box; 800/Case
181026IDB	1ml 26G x 3/8" IDB	100/Inner Box; 800/Case
181025S	1ml 25G x 5/8" Slip tip	100/Inner Box; 800/Case

1ml SOL-M™ Allergy Syringe Tray

TRAY PACKAGE

181026T	1ml 26G x 1/2"	25/Inner Box; 1000/Case
181126IDBT	1ml 26G x 1/2" IDB	25/Inner Box; 1000/Case
181026IDBT	1ml 26G x 3/8" IDB	25/Inner Box; 1000/Case
181027T	1ml 27G x 1/2"	25/Inner Box; 1000/Case
181027IDBT	1ml 27G x 3/8" IDB	25/Inner Box; 1000/Case
181027SNT	1ml 27G x 3/8"	25/Inner Box; 1000/Case

1ml SOL-M™ Allergy Mixing Tray

TRAY PACKAGE

181023T	1ml 23G x 1/2"	25/Inner Box; 1000/Case
---------	----------------	-------------------------

Please contact us for the availability in your country.



REF	Description	Qty
-----	-------------	-----




1ml SOL-M™ Syringe Only

BLISTER PACKAGE

180001	1ml	Luer Lock Tip	100/Inner Box; 800/Case
180011	1ml	Silp Tip	100/Inner Box; 800/Case

1ml SOL-M™ Luer Lock Syringe with Exchangeable Needle

BLISTER PACKAGE

18012912	 1ml	29G x 1/2"	100/Inner Box; 800/Case
012912PC	 1ml	29G x 1/2" (PC)	100/Inner Box; 800/Case
18012010	 1ml	20G x 1"	100/Inner Box; 800/Case













3ml SOL-M™ Syringe Only

BLISTER PACKAGE

180003	3ml	Luer Lock Tip	100/Inner Box; 800/Case
180003ST	3ml	Silp Tip	100/Inner Box; 800/Case

3ml SOL-M™ Luer Lock Syringe with Exchangeable Needle

BLISTER PACKAGE

1832558	 3ml	25G x 5/8"	100/Inner Box; 800/Case
1832510	 3ml	25G x 1"	100/Inner Box; 800/Case
1832515	 3ml	25G x 1 1/2"	100/Inner Box; 800/Case
1832310	 3ml	23G x 1"	100/Inner Box; 800/Case
1832210	 3ml	22G x 1"	100/Inner Box; 800/Case
18322125	 3ml	22G x 1 1/4"	100/Inner Box; 800/Case
1832215	 3ml	22G x 1 1/2"	100/Inner Box; 800/Case
1832110	 3ml	21G x 1"	100/Inner Box; 800/Case
1832115	 3ml	21G x 1 1/2"	100/Inner Box; 800/Case
1832010	 3ml	20G x 1"	100/Inner Box; 800/Case
1832015	 3ml	20G x 1 1/2"	100/Inner Box; 800/Case
1831815	 3ml	18G x 1 1/2"	100/Inner Box; 800/Case

5ml SOL-M™ Syringe Only

BLISTER PACKAGE

180005	5ml	Luer Lock Tip	100/Inner Box; 800/Case
180005ST	5ml	Silp Tip	100/Inner Box; 800/Case

REF	Description	Qty
-----	-------------	-----

5ml SOL-M™ Luer Lock Syringe with Exchangeable Needle

BLISTER PACKAGE

1852210	 5ml	22G x 1"	100/Inner Box; 800/Case
1852215	 5ml	22G x 1 1/2"	100/Inner Box; 800/Case
1852110	 5ml	21G x 1"	100/Inner Box; 800/Case
1852115	 5ml	21G x 1 1/2"	100/Inner Box; 800/Case
1852010	 5ml	20G x 1"	100/Inner Box; 800/Case
1852015	 5ml	20G x 1 1/2"	100/Inner Box; 800/Case




10ml SOL-M™ Syringe Only

BLISTER PACKAGE

180010	10ml	Luer Lock Tip	100/Inner Box; 1200/Case
180010ST	10ml	Silp Tip	100/Inner Box; 1200/Case
180010ET	10ml	Eccentric Tip	100/Inner Box; 1200/Case

10ml SOL-M™ Luer Lock Syringe with Exchangeable Needle

BLISTER PACKAGE

1812210	 10ml	22G x 1"	100/Inner Box; 1200/Case
1812110	 10ml	21G x 1"	100/Inner Box; 1200/Case
1812115	 10ml	21G x 1 1/2"	100/Inner Box; 1200/Case
1812010	 10ml	20G x 1"	100/Inner Box; 1200/Case
1812015	 10ml	20G x 1 1/2"	100/Inner Box; 1200/Case

20ml SOL-M™ Syringe Only

BLISTER PACKAGE

180020	20ml	Luer Lock Tip	100/Inner Box; 800/Case
180020ST	20ml	Luer Slip	100/Inner Box; 800/Case
180020ET	20ml	Eccentric Tip	100/Inner Box; 800/Case

30ml SOL-M™ Syringe Only

BLISTER PACKAGE

180030	30ml	Luer Lock Tip	50/Inner Box; 400/Case
180030ET	30ml	Eccentric Tip	50/Inner Box; 400/Case

60ml SOL-M™ Syringe Only

BLISTER PACKAGE

180060	60ml	Luer Lock Tip	30/Inner Box; 240/Case
180060ET	60ml	Eccentric Tip	30/Inner Box; 240/Case
180060CT	60ml	Catheter Tip	30/Inner Box; 240/Case

Please contact us for the availability in your country.

SOL-M™

Hypodermic Needle

SOL-Millennium offers a wide range of Hypodermic Needles.

Features & Benefits

Ultra sharp needle - Tricut bevel minimizes pain on injection

Translucent hubs has visual flashback confirmation

The hubs are color coded for gauge identification



REF	Description	Qty
SOL-M™ Hypodermic Needle		
113005	30 G x 1/2" regular bevel	100/Box; 1000/Case
112905	29 G x 1/2" regular bevel	100/Box; 1000/Case
112705	27 G x 1/2" regular bevel	100/Box; 1000/Case
112734	27 G x 3/4" regular bevel	100/Box; 1000/Case
112638	26 G x 3/8" IDB bevel	100/Box; 1000/Case
112612	26 G x 1/2" regular bevel	100/Box; 1000/Case
112658	26 G x 5/8" regular bevel	100/Box; 1000/Case
112506	25 G x 5/8" regular bevel	100/Box; 1000/Case
112510	25 G x 1" regular bevel	100/Box; 1000/Case
112515	25 G x 1 1/2" regular bevel	100/Box; 1000/Case
112410	24 G x 1" regular bevel	100/Box; 1000/Case
112334	23 G x 3/4" regular bevel	100/Box; 1000/Case
112310	23 G x 1" regular bevel	100/Box; 1000/Case
1123125	23 G x 1 1/4" regular bevel	100/Box; 1000/Case
112210	22 G x 1" regular bevel	100/Box; 1000/Case
112215	22 G x 1 1/2" regular bevel	100/Box; 1000/Case
1122125	22 G x 1 1/4" regular bevel	100/Box; 1000/Case
112110	21 G x 1" regular bevel	100/Box; 1000/Case
112115	21 G x 1 1/2" regular bevel	100/Box; 1000/Case
112010	20 G x 1" regular bevel	100/Box; 1000/Case
112015	20 G x 1 1/2" regular bevel	100/Box; 1000/Case
111915	19 G x 1 1/2" regular bevel	100/Box; 1000/Case
111810	18 G x 1" regular bevel	100/Box; 1000/Case
111815	18 G x 1 1/2" regular bevel	100/Box; 1000/Case

Please contact us for the availability in your country.



/ Diabetes




SOL-GUARD™ Insulin Syringe



SOL-GUARD™ Insulin Syringe, engineered safety mechanism to potentially reduce the risk of needlestick injury. The SOL-GUARD™ Insulin Syringe, has a transport mode, clear graduations for visibility of draw-up of medication and low dead space.

Features

- **Special Design Safety Shield**- Protection to healthcare worker during transporting device, and after use. Easy to lock either right or left twist, confirmation for final locking position is an audible click
- **Transport Position** - Potentially reducing the risk of needle stick injury while moving locations
- **Low Dead Space** - Minimizes medication waste
- **Clear, Bold Graduations** - Accuracy in draw up of medication



REF	Description	Qty
0.3mL SOL-GUARD™ Insulin Syringe		
BLISTER PACKAGE		
200003SG	0.3ml 29G x 1/2"	100/Box; 800/Case
0.5mL SOL-GUARD™ Insulin Syringe		
BLISTER PACKAGE		
200002SG	0.5ml 29G x 1/2"	100/Box; 800/Case
1mL SOL-GUARD™ Insulin Syringe		
BLISTER PACKAGE		
200017SG	1ml 29G x 1/2"	100/Box; 800/Case

Please contact us for the availability in your country.

SOL-CARE™

Insulin Safety Syringes with Fixed Needle



Sol-Millennium offers innovative products that bring comfort and convenience for insulin drug therapy. SOL-CARE™ Insulin Safety Syringes have low Dead Space on both our 100U and 50U.

Features & Benefits

- Special designed plunger with fit tip reduces dead space in syringe
- Ultra sharp needle minimizes pain on injection
- Patent Pending Locking Ring secures the needle into the barrel of the syringe

Easy to use


Fingers remain behind the needle during activation of safety mechanism

Syringe cannot be reused

Clear, bold gradations

Visual and audible confirmation of safety activation

Potential reduction in sharps disposal costs



REF	Description	Qty
0.5ml SOL-CARE™ Insulin Safety Syringe with Fixed Needle		
BLISTER PACKAGE		
1000021M	0.5ml 29G x 1/2"	100/Box; 1000/Case
1000891M	0.5ml 30G x 5/16"	100/Box; 1000/Case
1ml SOL-CARE™ Insulin Safety Syringe with Fixed Needle		
BLISTER PACKAGE		
1000711M	1ml 28G x 1/2"	100/Box; 1000/Case
1000171M	1ml 29G x 1/2"	100/Box; 1000/Case
1000201M	1ml 30G x 5/16"	100/Box; 1000/Case
1000811M	1ml 30G x 1/2"	100/Box; 1000/Case

Please contact us for the availability in your country.

SOL-M™


Insulin Syringe with Fixed Needle

Sol-Millennium offers innovative products that bring comfort and convenience for diabetes drug therapy. The package is tamper evident for protection of the sterile contents.

Features

Our full line of SOL-M™ Insulin Syringe adapted to the individual needs and preferences of insulin needs. A wide range of needle sizes.

- **Low dead space** - Minimizes residual medication wastage.
- **Optimal needle sharpness and smooth plunger movement** - Makes injections more comfortable.
- **Transparent barrel with bold scale** - Markings make it easier to see and measure the syringe content.
- **Flat top rubber stopper** - For easy lining up with scale marks.



REF	Description	Qty
-----	-------------	-----

0.3ml SOL-M™ Insulin Syringe with Fixed Needle

SELF-CONTAINED(10pcs in 1 PE bag)

1633151B	0.3ml 31G x 5/16"	100/Inner Box; 1000/Case
1633051B	0.3ml 30G x 5/16"	100/Inner Box; 1000/Case
1633012B	0.3ml 30G x 1/2"	100/Inner Box; 1000/Case
1632912B	0.3ml 29G x 1/2"	100/Inner Box; 1000/Case

BLISTER PACKAGE

16330516	0.3ml 30G x 5/16"	100/Inner Box; 800/Case
1633012	0.3ml 30G x 1/2"	100/Inner Box; 800/Case
1632912	0.3ml 29G x 1/2"	100/Inner Box; 800/Case

0.5ml SOL-M™ Insulin Syringe with Fixed Needle

SELF-CONTAINED(10pcs in 1 PE bag)

1653151B	0.5ml 31G x 5/16"	100/Inner Box; 1000/Case
1653051B	0.5ml 30G x 5/16"	100/Inner Box; 1000/Case
1653012B	0.5ml 30G x 1/2"	100/Inner Box; 1000/Case
1652912B	0.5ml 29G x 1/2"	100/Inner Box; 1000/Case

BLISTER PACKAGE

16530516	0.5ml 30G x 5/16"	100/Inner Box; 800/Case
1653012	0.5ml 30G x 1/2"	100/Inner Box; 800/Case
1652912	0.5ml 29G x 1/2"	100/Inner Box; 800/Case

1ml SOL-M™ Insulin Syringe with Fixed Needle

SELF-CONTAINED(10pcs in 1 PE bag)

1613151B	1ml 31G x 5/16"	100/Inner Box; 1000/Case
1613051B	1ml 30G x 5/16"	100/Inner Box; 1000/Case
1613012B	1ml 30G x 1/2"	100/Inner Box; 1000/Case
1612912B	1ml 29G x 1/2"	100/Inner Box; 1000/Case

BLISTER PACKAGE

16130516	1ml 30G x 5/16"	100/Inner Box; 800/Case
1613012	1ml 30G x 1/2"	100/Inner Box; 800/Case
1612912	1ml 29G x 1/2"	100/Inner Box; 800/Case

Please contact us for the availability in your country.

SOL-M™


Pen Needle



SOL-M™ Pen Needle is your simple and natural choice for daily diabetes drug therapy. They can be universally applied to most insulin pens on the market providing great convenience for your home injection. SOL-M™ Pen Needles are available in a variety of sizes, offering you perfect options for your own requirement.

Features

- **Premium needle quality and extra needle sharpness** - Minimizes patient discomfort.
- **Screw thread** - Firm and secure fit for the cannula.
- **Thin wall technology** - Ensures optimal insulin flow.
- **New 4mm 32G pen needle** - Effective delivery of insulin & greater injection comfort.



REF	Description	Qty
SOL-M™ Pen Needle		
BLISTER PACKAGE		
N42006	29G x 12.5mm	100/Inner Box; 1200/Case
N42003	31G x 8mm	100/Inner Box; 1200/Case
N42004	31G x 5mm	100/Inner Box; 1200/Case
N42005	31G x 4mm	100/Inner Box; 1200/Case
N42002	31G x 6mm	100/Inner Box; 1200/Case
N42001	30G x 8mm	100/Inner Box; 1200/Case
N42007	32G x 4mm	100/Inner Box; 1200/Case

Please contact us for the availability in your country.



/ Allergy



SOL-CARE™

Safety Syringe Allergy Tray

SOL-CARE™ Safety Syringes are designed to reduce the risk of needlestick injury to healthcare workers. The needle is retracted into the barrel of the syringe and secured with our patent pending locking ring eliminating syringe reuse.



Features & Benefits

- Special designed plunger with fit tip reduces dead space in syringe
- Ultra sharp needle minimizes pain on injection
- Patent Pending Locking Ring secures the needle into the barrel of the syringe
- Easy to use
- Fingers remain behind the needle during activation of safety mechanism
- Syringe cannot be reused
- Clear, bold gradations
- Visual and audible confirmation of safety activation
- Potential reduction in sharps disposal costs

REF	Description	Qty
1ml SOL-CARE™ Safety Syringe Allergy Tray		
TRAY PACKAGE		
100032IM	■ 1ml 27G x 1/2"	25/Tray; 1000/Case
100069IM	■ 1ml 26G x 3/8" IDB	25/Tray; 1000/Case

Please contact us for the availability in your country.



/ Immunization



SOL-CARE™

Safety Syringe Flu Tray

SOL-CARE™ Safety Syringes are designed to reduce the risk of needlestick injury to healthcare workers. The needle is retracted into the barrel of the syringe and secured with our patent pending locking ring eliminating syringe reuse.



Features & Benefits

- S**pecial designed plunger with fit tip reduces dead space in syringe
- U**ltra sharp needle minimizes pain on injection
- P**atent Pending Locking Ring secures the needle into the barrel of the syringe
- E**asy to use
- F**ingers remain behind the needle during activation of safety mechanism
- S**yringe cannot be reused
- C**lear, bold gradations
- V**isual and audible confirmation of safety activation
- P**otential reduction in sharps disposal costs

REF	Description	Qty
1ml SOL-CARE™ Safety Syringe Flu Tray		
TRAY PACKAGE		
100088IM	1ml 25G x 1"	25/Tray; 1000/Case

Please contact us for the availability in your country.

SOL-CARE™

Safety Syringe for Immunization

An evolutionary safety solution designed for improving the efficiency of Immunization injection.



Features & Benefits

- S**pecial designed plunger with fit tip reduces dead space in syringe
- U**ltra sharp needle minimizes pain on injection
- P**atent Pending Locking Ring secures the needle into the barrel of the syringe
- E**asy to use
- F**ingers remain behind the needle during activation of safety mechanism
- S**yringe cannot be reused
- C**lear, bold gradations
- V**isual and audible confirmation of safety activation
- P**otential reduction in sharps disposal costs

REF	Description	Qty
0.5ml SOL-CARE™ Safety Syringe with Fixed Needle		
BLISTER PACKAGE		
100007IM	0.5ml 25G x 1"	100/Box; 1000/Case
1ml SOL-CARE™ Safety Syringe with Fixed Needle		
BLISTER PACKAGE		
100018IM	1ml 25G x 5/8"	100/Box; 1000/Case
100034IM	1ml 25G x 1"	100/Box; 1000/Case

Please contact us for the availability in your country.

SOL-Millennium® Safety Box

Safety Box Is for safe and easy handling and disposal of used syringes and needles. Proper use of Safety Box can help prevent virus transmission by needlestick injuries.

EASY

- * Pictorial instructions and squeeze code on flaps for easy assembling.
- * Clear visibility of maximum fill line.
- * Easy to carry for outside vaccination.
- * Combustible without lighter fluid.

SAFE

- * Handle embodied to the upper side allows for one hand handling.
- * Integrated sealing strip prevent spillage of contents.
- * Prevents health workers and communities from being exposed to contaminated needles.
- * On site destruction prevent downstream reuse of disposable syringes.
- * Waterproof and resistant to punctures.

EFFICIENCY

- * Can be used anywhere (Vaccination campaigns, Hospitals, Mobile health clinics, Medic in the armed forces).

SPECIFICATION

Volume :	5 Liters
External dimension	
Before assembling :	590 x 300 x 4mm
After assembling :	280 x 162 x 122mm
Material :	recycled board, waterproof, and puncture proof
Printing :	brown background with black illustrations and texts
Thickness of wall :	1.2mm
Weight fully assembled :	280grams
Diam. syringe insert hole :	38mm



REF	Description	Qty
SOL-Millennium™ Safety Box		
SX30500	5L Safety Box	25/Case

Please contact us for the availability in your country.

A close-up photograph of a hand holding a small glass vial and another hand holding a syringe. The syringe is being used to draw a yellowish liquid from the vial. The background is a white grid pattern with green and pink curved shapes on the left side.

/ Pharmacy


SOL-M™ Blunt Fill Needle & SOL-M™ Blunt Fill Needle with 5 micron filter

SOL-M™ Blunt Fill Needle is designed to potentially reduce the risk of needlestick injury while drawing up medication, medication re-constitution, and accessing IV Bags. SOL-M™ Blunt Fill Needles address a wide range of medication preparation needs. The SOL-M™ Blunt Fill Needle provides a substitute for a conventional needle reducing the risk of needle stick injuries.

SOL-M™ Blunt Fill Needle with filter features a 5 micron filter membrane that is designed to prevent glass particles being drawn into the syringe when preparing medications from ampules.

Features & Benefits

- 45° bevel reducing the risk of needle stick injury
 - 18G for ease in draw-up of medications
 - Traditional hub that fits luer lock or slip tip syringes
 - 5 micron filter—Effective filtration of glass particles
- Color coded for identification
Lower cost than sharps engineered devices



REF	Description	Qty
SOL-M™ Blunt Fill Needle		
110022	■ 18 x 1 1/2"	100/Box; 1000/Case
110022F	■ 18 x 1 1/2" with 5 micron filter	100/Box; 1000/Case
110023*	■ 18 x 50mm	100/Box; 1000/Case
110023F*	■ 18 x 50mm with 5 micron filter	100/Box; 1000/Case
* Not available in US		

Please contact us for the availability in your country.

A close-up photograph of a person's hands wearing white latex gloves. The hands are holding a clear plastic syringe. The right hand is at the top, holding a small purple cap that is being removed from the top of the syringe. The left hand is at the bottom, holding the main body of the syringe. The syringe has a black plunger and a clear barrel with markings. The background is white. On the left side of the image, there are decorative curved shapes in shades of purple and blue. On the right side, there is a solid purple shape.

/ **S**afety Products for Nuclear Medicine

SOL-CARE™


Luer Lock Safety Syringe

—Fits most Nuclear Medicine Shields

Our 3ml Luer Lock Safety Syringe is designed to be used and fits most Nuclear Medicine Shields for protection during draw-up and administration of medication. A purple protective cap is placed on the open end of the syringe which provides additional safety after completion of the injection and activation of the safety mechanism.

Features & Benefits

- All SOL-CARE™ safety syringes are now enhanced with our patented no gap design allowing clinicians to push plunger to zero line when injecting air into a vial
- Exchangeable needles allow clinician to change needle after draw-up of medication
- Special designed plunger with fit tip reduces dead space in syringe
- Patent Pending Locking Ring secures the needle in the barrel of the syringe preventing re-use of syringe
- Compatible with all standard hypodermic needles
- Easy to use
- Fingers remain behind the needle during activation of safety mechanism
- Clear, bold gradations
- Visual and audible confirmation of safety activation
- Potential reduction in sharps disposal costs
- Ultra sharp tricut bevel needle minimizes pain on injection



REF	Description	Qty
SOL-CARE™ Luer Lock Safety Syringe – Fits in most nuclear medicine shields		
BLISTER PACKAGE		
150085IM	■ 3ml 22G x 1"	100/Box; 800/Case
150085IM-C	Purple Protective Cap	100/Box; 800/Case

Please contact us for the availability in your country.

/ **B**lood Collection



SOL-CARE™

Safety Blood Collection Tube Holder Set



Features & Benefits

- * The SOL-CARE™ Safety Blood Collection Tube Holder Set features a Blood Collection Tube Holder and a Vacuum Activation Tube. The safety mechanism is an integral part of the device.
- * The activation virtually eliminates exposure to both ends of the contaminated needle, reducing the risk of needlestick injury.
- * The SOL-CARE™ Blood Collection Tube Holder is compatible with standard, multi-sample blood collection needles.
- * The SOL-CARE™ Safety Blood Collection Tube Holder Set allows users to keep both hands behind the needle when activating the retraction of the needle.
- * The Single Use Safety Blood Collection Tube Holder Set complies with OSHA regulations to help eliminate cross contamination.



REF	Description	Qty
SOL-CARE™ Safety Blood Collection Tube Holder Set		
BLISTER PACKAGE		
7001BCD	Safety Blood Collection Tube Holder Set	25/Box; 200/Case

Please contact us for the availability in your country.



/ Oral Feeding

SOL-M™

Oral Dispensing Syringe



Sol-millennium offers a broad size of oral dispensing syringes including both clear and amber syringes.

Features

- **Special tip design** - for wrong-route administration prevention
- **Clear dose scale** - for accurate volume confirmation
- **Amber barrel** - design for protect light-sensitive drugs

REF	Description	Qty
SOL-M™ Oral Dispensing Syringe		
Clear with Tip Cap		
21001OC	1 mL Oral Syringe	100/Inner Box; 400/Case
21003OC	3 mL Oral Syringe	100/Inner Box; 400/Case
21005OC	5 mL Oral Syringe	100/Inner Box; 400/Case
21010OC	10 mL Oral Syringe	100/Inner Box; 400/Case
Clear with Tip Cap with Purple Plunger		
21001OP	1 mL Oral Syringe	100/Inner Box; 400/Case
21003OP	3 mL Oral Syringe	100/Inner Box; 400/Case
21005OP	5 mL Oral Syringe	100/Inner Box; 400/Case
21010OP	10 mL Oral Syringe	100/Inner Box; 400/Case
Amber with Tip Cap		
21001OA	1 mL Oral Syringe	100/Inner Box; 400/Case
21003OA	3 mL Oral Syringe	100/Inner Box; 400/Case
21005OA	5 mL Oral Syringe	100/Inner Box; 400/Case
20010OA	10 mL Oral Syringe	100/Inner Box; 400/Case

SOL-M™ Oral Dispensing Syringe

Clear with Tip Cap

21001OC	1 mL Oral Syringe	100/Inner Box; 400/Case
21003OC	3 mL Oral Syringe	100/Inner Box; 400/Case
21005OC	5 mL Oral Syringe	100/Inner Box; 400/Case
21010OC	10 mL Oral Syringe	100/Inner Box; 400/Case


Clear with Tip Cap with Purple Plunger

21001OP	1 mL Oral Syringe	100/Inner Box; 400/Case
21003OP	3 mL Oral Syringe	100/Inner Box; 400/Case
21005OP	5 mL Oral Syringe	100/Inner Box; 400/Case
21010OP	10 mL Oral Syringe	100/Inner Box; 400/Case

Amber with Tip Cap

21001OA	1 mL Oral Syringe	100/Inner Box; 400/Case
21003OA	3 mL Oral Syringe	100/Inner Box; 400/Case
21005OA	5 mL Oral Syringe	100/Inner Box; 400/Case
20010OA	10 mL Oral Syringe	100/Inner Box; 400/Case

Please contact us for the availability in your country.






/ **P**et Care



SOL-M™

Insulin Syringe with Fixed Needle For U40 and U100 Insulin




Sol-Millennium offers traditional diabetic syringes for the pet care market.

Features & Benefits

SOL-M™ Traditional Insulin Syringe is specifically designed for pets

- **Low dead space** - Minimizes residual medication wastage
- **Optimal needle sharpness and smooth plunger movement** - Makes injections more comfortable
- **Transparent barrel with bold scale** - Markings make it easier to see and measure the syringe content
- **Flat top rubber stopper** - For more accurate dosing



REF	Description	Qty
0.3ml SOL-M™ U40 Insulin Syringe with Fixed Needle		
SELF-CONTAINED(10pcs in 1 PE bag)		
V32905-40	0.3ml 29G x 1/2"	100/Box; 1000/Case
V3305-40	0.3ml 30G x 1/2"	100/Box; 1000/Case
0.5ml SOL-M™ U40 Insulin Syringe with Fixed Needle		
SELF-CONTAINED(10pcs in 1 PE bag)		
V5295-40	0.5ml 29G x 1/2"	100/Box; 1000/Case
V5305-40	0.5ml 30G x 1/2"	100/Box; 1000/Case
1ml SOL-M™ U40 Insulin Syringe with Fixed Needle		
SELF-CONTAINED(10pcs in 1 PE bag)		
V12905-40	1ml 29G x 1/2"	100/Box; 1000/Case
V1305-40	1ml 30G x 1/2"	100/Box; 1000/Case

REF	Description	Qty
0.3ml SOL-M™ U100 Insulin Syringe with Fixed Needle		
SELF-CONTAINED(10pcs in 1 PE bag)		
V3295	0.3ml 29G x 1/2"	100/Box; 1000/Case
V3305	0.3ml 30G x 1/2"	100/Box; 1000/Case
0.5ml SOL-M™ U100 Insulin Syringe with Fixed Needle		
SELF-CONTAINED(10pcs in 1 PE bag)		
V5295	0.5ml 29G x 1/2"	100/Box; 1000/Case
V5305	0.5ml 30G x 1/2"	100/Box; 1000/Case
1ml SOL-M™ U100 Insulin Syringe with Fixed Needle		
SELF-CONTAINED(10pcs in 1 PE bag)		
V12905	1ml 29G x 1/2"	100/Box; 1000/Case
V1305	1ml 30G x 1/2"	100/Box; 1000/Case

Please contact us for the availability in your country.

SOL-MILLENNIUM MEDICAL

Shanghai Office:

Room 2202, No. 58 Changliu Road
Pudong, Shanghai, China, 200135
Tel: (86 21) 6862 2586
Fax: (86 21) 6862 2585
E-mail: info@sol-m.com
Web: www.sol-m.com

North America Office:

5415 Sugarloaf Parkway, Suite 2203
Lawrenceville, GA 30043
Tel: 404.973.2200
Fax: 404.973.2212
E-mail: info@sol-m.com
Web: www.sol-millenniumna.com

Europe Office:

Emilii Plater 53, PL-00113 Warszawa
Tel: +48 22 528 69 41
Fax: +48 22 528 67 01
E-mail: info@sol-m.com
Web: www.sol-millenniumeu.com

Australia Office:

Suite 404/8-12 Station St
Homebush, NSW, 2140
Tel: 02 9746 7738
E-mail: info@sol-m.com
Web: www.sol-m.com

Brasil Office:

Avenida Anhanguera , 869
Ribeirão Preto - São Paulo, Brasil
Tel: 14.025-480
Telefone: 55 16 3442 1454 / 55 11 984060126
E-mail: info@sol-m.com
Web: www.sol-m.com